

THE TOP



ACCURATE CESIUM CHLORIDE

IT WON'T SPIN YOU BROKE!

A1059	CESIUM CHLORIDE PREP™	1 kg	\$127.50
A1250	CESIUM CHLORIDE SPECIAL BIO.™	100 gm	37.50
		1 kg	177.00
A1255	CESIUM CHLORIDE OPTI-MAX.™	100 gm	52.50
	0.01 MAX. DENSITY GRADIENT	1 kg	350.00

F.O.B. WESTBURY, NEW YORK TERMS: 2% 10 NET 30 DAYS

 **ACCURATE Chemical & Scientific Corporation**

300 SHAMES DRIVE, WESTBURY, NY 11590 • TEL: 516-433-4900 • TELEX: 14-4617


WEST COAST SALES:

714-235-9400

NEW YORK CITY:

212-724-7800

TO ORDER CALL TOLL FREE **1-800-645-6264** (Except N.Y. State call 516-433-4900)



***Serums—
from 15 years
of tissue culture experience.***

***Quality — from a century of
pharmaceutical discipline.***

Armour Pharmaceutical has been a pharmaceutically oriented group since our founding in 1886. Because of this orientation, we have the "pharmaceutical attitude" generated over the years by our self-imposed quality standards and those associated with changing regulatory demands.

We are more than just a biochemical manufacturer — we are a pharmaceutical company.

*Tissue Culture Reagents
REHATUIN®
NORMAL CALF SERUM
HORSE SERUM*

*Sterile Trypsin Solution
ENZAR® T*



Armour Pharmaceutical Company

Biochemicals — a commitment to the '80s



Armour Pharmaceutical Tissue Culture Reagents are of the highest quality available! Our production, R & D and QA facilities in Kankakee, Illinois are absolutely "state-of-the-science" modern. Plus, we retain the most qualified consultants to help keep us totally up-to-date on both processing and testing methods.

We are interested in answering your needs for serums. We pre-test to be sure that what you get is consistent with what you want.

In addition, we offer complete technical service as a valued and trustworthy resource to our customers.

REHATUIN®

(Fetal Bovine Serum Sterile Filtered)
Rehatuin is a sterile filtered fetal bovine serum collected by aseptic methods.

This aseptic collection minimizes the chances of bacterial contamination and resultant microbial endotoxins. Fetal bovine blood is collected in federally inspected packing houses. To insure the highest quality, the blood is immediately cooled, clotted and centrifuged. The pooled serum is pretested for growth promotion using at least four different cell lines and screened for the presence of mycoplasma and adventitious viral agents prior to sterile filtration. Rehatuin undergoes the same stringent USP XX sterility testing used in pharmaceutical operations. Rehatuin also meets federal regulations (9CFR-113.53) and is pre-tested to meet the needs of your unique cell lines.

The quality of your research is directly related to, the quality of the reagent you use.

NORMAL CALF SERUM

Normal calf serum is a sterile filtered serum derived from specially fed calves maintained on a special diet according to our specifications designed to yield serum with optimal growth properties. The serum is collected in federally inspected packing houses. The collection and rapid processing minimize the chances of bacterial contamination and resultant microbial endotoxins.

Normal calf serum is a low cost additive to tissue culture media for cell growth support. Our normal calf serum has been preselected to yield growth parameters equivalent to those obtained with fetal bovine serum in at least four different cell lines. We routinely pre-test the serum to meet your unique cell line needs.

HORSE SERUM

Horse serum is collected from horses approved by USDA veterinarians who assure the animals are healthy and free of infectious equine anemia, glanders and piroplasmiasis.

The blood is collected into sterile

containers and immediately chilled and centrifuged. The serum is sterile filtered through a final filter and packaged.

Horse serum is used in tissue culture media for growth promotion and is pre-tested to meet your unique cell line needs. Horse serum is also a valuable additive in the propagation of mycoplasma and bacteria in microbiology laboratories.

ENZAR®-T

(sterile Trypsin solution)

ENZAR®-T is a Trypsinizing concentrate prepared from porcine pancreatic enzymes. It contains enzymes naturally occurring in porcine pancreas (trypsin, chymotrypsin and elastase) and offers product uniformity never before attained.

ENZAR®-T solution is treated by a process which has been demonstrated to inactivate mycoplasma, bacteria, fungi and viruses. It retains its utility at normal refrigerator temperatures for at least three to four months.

ENZAR®-T solution can be used with excellent results in any cell system requiring Trypsinization.

For an interchange of technical information on Tissue Culture Reagents and Trypsin solutions, call 1-800-435-1852 (in Illinois 1-312-726-6851) or write: Technical Service Manager, Biochemicals Armour Pharmaceutical Company, P.O. Box 511, Kankakee, IL 60901



Armour Pharmaceutical Company

303 South Broadway, Tarrytown, NY 10591 • 1-914-631-8888

Biochemicals – a commitment to the '80s

The newest addition to the annual
MICROBIOLOGY series,

MICROBIOLOGY-

1982

*puts into focus
exciting new work in
the following major areas of microbiological research:*

CELL-CELL INTERACTIONS IN PATHOGENESIS

Bacterial Adhesion in Pathogenesis

- Introduction
- Autochthonous Populations
- Role of Surface Proteins Involved in Surface Adhesion
- Role of Surface Polysaccharides in Bacterial Adhesion
- Variable Surface Factors in Bacterial Adhesion

The Macrophage in Host Defense

Cell Surface Receptors Involved in the Immune System

COMPARATIVE MOLECULAR BIOLOGY

Streptococcal Genetics

- Conjugation
- Transformation, Transduction, and Bacteriophage
- Drug Resistance
- Genetics and Physiology
- Metabolic and Bacteriocin-Related Plasmids
- Development and Use of Recombinant DNA Technology

Bacillus subtilis Chromosome

- Structure and Modification
- Molecular Cloning

As usual in this well-established series, these extended abstracts are carefully selected for their importance to the development of the knowledge base in each of these areas. The book will be particularly useful in research laboratories and institutional libraries.

Bibliographic and Ordering Information

Editor: David Schlessinger.

July 1982, 400 pages, clothbound. \$15.00 to ASM members; \$25.00 to others. Please note that ASM members must remit payment with order. Orders from outside the United States and Canada must be prepaid before shipment. Prices are subject to change without notice. ISBN: 0-914826-42-5. Standing orders to the series are welcome.

Previous volumes in the series, which started publication in 1975 with the volume entitled *Microbiology — 1974*, are still available in limited quantities. Prices for each preceding volume are \$15.00 for ASM members and \$25.00 for others.

ASM
American Society
for Microbiology

Published and distributed by:
American Society for Microbiology
Publications Office
1913 I Street, NW
Washington, DC 20006

Meeting the bacteriologist's need for a compact handbook of laboratory methods—

Manual of Methods for General Bacteriology

EDITOR-IN-CHIEF:
Philipp Gerhardt

SECTION EDITORS:
Ralph N. Costilow
Noel R. Krieg
R. G. E. Murray
Eugene W. Nester
G. Briggs Phillips
Willis A. Wood

Published and distributed by:



American Society for Microbiology
1913 I Street, NW
Washington, DC 20006

This *Manual* covers all the taxa of bacteria, and its approach is generalized for both pure and applied bacteriology. Aimed at senior undergraduate and beginning graduate science majors, at teachers and researchers, the *Manual* provides practical descriptions of laboratory procedures that are known to be reliable and that are in frequent use. It is an invaluable tool for all involved in bacteriology in the laboratory setting.

The book comprises six sections, each of which is self-contained, with its own table of contents to help the reader find specific methods as easily as possible. The contents are listed below.

METHODOLOGY FOR GENERAL BACTERIOLOGY

MORPHOLOGY—Light Microscopy—Specimen Preparation for Light Microscopy—Determinative Methods of Light Microscopy—Electron Microscopy—Cell Fractionation

GROWTH—Biophysical Factors in Growth—Biochemical Factors in Growth—Enrichment and Isolation—Solid Culture—Liquid Culture—Growth Measurement—Preservation

GENETICS—Gene Mutation—Gene Transfer—Plasmids

METABOLISM—Physical Methods—Chemical Composition—Enzymatic Activity—Permeability and Transport

SYSTEMATICS—General Characterization—Numerical Taxonomy—Genetic Characterization

LABORATORY SAFETY—Sterilization—Containment and Disinfection

ADDENDUM—Diluents and Biomass Measurement

ORDERING INFORMATION

Publication date: March 1981. 525 pages, fully illustrated.

Clothbound


Flexible binding:

ISBN: 0-914826-29-8

ISBN: 0-914826-30-1

\$28.00 (\$20.00 to members of ASM) \$24.00 (\$17.00 to members of ASM)

Prices subject to change without notice.



I am a member of ASM.

Send me _____ copy(ies)

(limit of 3 at member prices)

\$20.00 (Cloth);

\$17.00 (Flexible).

I have enclosed \$ _____.

Please send me _____ copy(ies) of

MANUAL OF METHODS FOR GENERAL BACTERIOLOGY

\$28.00 (Cloth);

\$24.00 (Flexible).

Send book(s) to:

I have enclosed \$ _____.

Name _____

Address _____

City/State/Zip _____

The successful new journal from ASM

Molecular and Cellular Biology

Editor-in-Chief: Aaron J. Shatkin

Editors: Harvey F. Lodish, Keith R. Porter, Louis Siminovitch, Paul S. Sypherd

Editorial Board: Bruce Adams, Guenter Albrecht-Buehler, Renato Baserga, Blair Bowers, Stuart Brody, Breck Byers, Lawrence Chasin, James Darnell, Rowland Davis, Gordon Dixon, Argiris Efstratiadis, Gary Felsenfeld, Werner Franke, James Friesen, Samuel Horowitz, Elias Lazarides, John B. Little, William Loomis, P.T. Magee, Harvey Ozer, Leo Parks, Mark Pearson, Jeremy Pickett-Heaps, Robert Pollack, John Pringle, Daniel B. Rifkin, Robert Roeder, David Sabatini, Phillip Sharp, Fred Sherman, Pamela Stanley, Joan Steitz, Walther Stoeckenius, James Van Etten, Bernard Weinstein, Harold Weintraub, Leslie Wilson

Scope Statement:

Molecular and Cellular Biology is devoted to the advancement and dissemination of fundamental knowledge concerning the molecular biology of eucaryotic cells, of both microbial and higher organisms. Papers on cellular morphology and function, genome organization, the regulation of genetic expression, morphogenesis, and somatic cell genetics are invited for submission.

Manuscripts:

Please submit manuscripts to the ASM Publications Office, 1913 I Street, NW, Washington, DC 20006.

Subscriptions:

Members may select MCB from the list of journal options on the Dues Bill. *Non-members* should send their order and \$70.00 for 12 issues to: M. Poe, American Society for Microbiology, PO Box 1192, Birmingham, AL 35201.

Published monthly.
About 1,600 pages a year.



Published by:
American Society for Microbiology
1913 I Street, NW
Washington, DC 20006

State-of-the-art technology in clinical microbiology

RAPID METHODS and AUTOMATION in Microbiology

EDITOR: Richard C. Tilton

This new collection of research papers on a topic of increasing and critical interest promises to be of real assistance to clinical microbiologists, clinical pathologists, immunologists, medical technologists, and others who are caught up in the rapidly changing environment of microbiology. The speed of publication lends added utility to this survey of technologically advanced—and advancing—methods and procedures. The book, indeed, captures some of the excitement that was in the air at the meeting in May 1981 on which this volume is based: the Third International Symposium on Rapid Methods and Automation in Microbiology.

As the keynote seminar of the symposium suggested, topics of particular concern include systems approaches, evaluation, and the future of diagnostic systems in microbiology. The book fully develops the following areas:

- Computer Systems in Clinical Microbiology
- Rapid Detection of Microorganisms in Body Fluids
- Hydrophobic Grid-Membrane Microbiology
- Rapid Reporting in Microbiology
- Immunological Diagnosis of Infectious Diseases
- Chromatographic Analysis: a New Future for Clinical Microbiology
- Rapid Detection of Microorganisms in Urine
- Rapid Methods in Food Microbiology
- Systems Approaches to Diagnostic Microbiology
- Rapid Identification of Viral, Chlamydial, and Mycoplasmal Agents
- Rapid Identification of Clinically Important Anaerobes
- Applications of the Autobac System
- Clinical and Research Applications of the MS-2 System
- Nonautomated Rapid Methods in Microbiology

ORDERING INFORMATION:

February 1982. 361 pages. Clothbound. \$30.00. ISBN: 0-914826-38-7

Price subject to change without notice.



Order from the publisher,

AMERICAN SOCIETY for MICROBIOLOGY

Publication Sales, 1913 I Street, NW / Washington, DC 20006

Better by Design

Contur™ Tissue Culture Dish



© 1982 Lab-Tek Division, Miles Laboratories, Inc.

Highlights

Consistent reproducible plating efficiency is assured through our well regulated manufacturing process. Technologically-advanced Lux Contur Dishes are made of virgin polystyrene for flat, distortion-free surfaces with exceptional optical clarity. Also available in a special plastic material, Permanox, for use in fixation and embedding for EM, where resistance to acetone and other reagents is required. Permanox also permits high transmission of UV light for studying light-influenced cells and in fluorescent studies where the excitation wave length is in the UV range.

Lux comprehensive line of Tissue Culture Ware includes Contur Dishes, standard style dishes and a Suspension Dish for use with non-adherent cell populations; as well as flasks, tubes, multiwell plates, and coverslips.

For information contact: Lab-Tek Division, Miles Laboratories, Inc., 30W475 North Aurora Road, Naperville, IL 60566. To order Lux® Tissue Culture Ware, contact our distributor: Flow Laboratories, McLean, Virginia.

Advancing
Technology
Committed to the
Tissue Culturist

Lab-Tek
Division

