



# **ACCURATE** **PNEUMOCOCCAL** **ANTISERA**

- **Omni-Serum**
- **Nine Pooled Sera**
- **46 Type or Group Sera**

---

**STEROID HORMONE**  
**ANTISERA**

---

**RADIOIMMUNOASSAY**  
**STANDARDS**

SEND FOR FREE LITERATURE

**ACCURATE** Chemical & Scientific Corporation

300 SHAMES DRIVE, WESTBURY, NY 11590 • TEL: 516-433-4900 • TELEX: 14-4817

WEST COAST SALES:  
619-235-9400

NEW YORK CITY:  
212-724-7800

TO ORDER CALL TOLL FREE 1-800-645-6264 (Except N.Y. State call 516-433-4900)

The newest addition to the annual  
MICROBIOLOGY series,

# **MICROBIOLOGY-**

# **1982**

puts into focus  
exciting new work in  
the following major areas of microbiological research:

## **CELL-CELL INTERACTIONS IN PATHOGENESIS**

Bacterial Adhesion in Pathogenesis

- Introduction
- Autochthonous Populations
- Role of Surface Proteins Involved in Surface Adhesion
- Role of Surface Polysaccharides in Bacterial Adhesion
- Variable Surface Factors in Bacterial Adhesion

The Macrophage in Host Defense

Cell Surface Receptors Involved in the Immune System

## **COMPARATIVE MOLECULAR BIOLOGY**

Streptococcal Genetics

- Conjugation
- Transformation, Transduction, and Bacteriophage
- Drug Resistance
- Genetics and Physiology
- Metabolic and Bacteriocin-Related Plasmids
- Development and Use of Recombinant DNA Technology

*Bacillus subtilis* Chromosome

- Structure and Modification
- Molecular Cloning

As usual in this well-established series, these extended abstracts are carefully selected for their importance to the development of the knowledge base in each of these areas. The book will be particularly useful in research laboratories and institutional libraries.

### **Bibliographic and Ordering Information**

Editor: David Schlessinger.

July 1982, 420 pages, clothbound. \$15.00 to ASM members; \$25.00 to others. Please note that ASM members must remit payment with order. Orders from outside the United States and Canada must be prepaid before shipment. Prices are subject to change without notice. ISBN: 0-914826-42-5. Standing orders to the series are welcome.

Previous volumes in the series, which started publication in 1975 with the volume entitled *Microbiology — 1974*, are still available in limited quantities. Prices for each preceding volume are \$15.00 for ASM members and \$25.00 for others.

**ASM**  
American Society  
for Microbiology

Published and distributed by:

**American Society for Microbiology**  
Publications Office  
1913 I Street, NW  
Washington, DC 20006

# NO EXCUSE FOR IGNORING SAFETY WHEN DNA SEQUENCING ON LONG GELS

**Savant's New High Voltage Package for every research laboratory that operates long gel cells in an exposed manner.**

**Model HV5000-150** — 5000 volt 150ma Power Supply including a 24 volt interlock electrical circuit \$1,475.00

**Model TTE32** — Table Top Enclosure Fire retardant, plastic enclosure with interlocked front and top doors 1,495.00

Original Price	\$2,970.00
<b>Safety Premium</b>	<b>500.00</b>
Total Price*	\$2,470.00


NOTE: Cables must be ordered separately see brochure.

\*Total "Safety Premium" Price available only when Power Supply and Enclosure are ordered at the same time.




## Safety Package Features:

- Operate 2 or more gel cells simultaneously
- Fire retardant plastic construction throughout
- Operate safely at 5KV
- 24 volt safety interlocked



## Improve Your Assays in Any Biochemical Procedure!

The Savant "NO BUMP" Speed Vac Concentrator/Evaporator is now standard procedure for the following and many other applications.



### Achieve 100% Sample Recovery!

- HPLC sample preparation.
- GLC sample preparation.
- DNA sequencing — evaporation of piperidine.
- PHOSPHOLIPID separation by HPLC — total evaporation of Acetonitrile-Water Solvent.
- TRIETHANOL AMINE carbonate evaporation.
- SOMATOSTATIN concentration by evaporation of acetic acid.
- INTERFERON purification by dialysis and evaporation of solvent phase.
- TRI-CYCLIC anti-depressant drug blood levels.
- NEUROTRANSMITTERS such as Acetylcholine, Dopa.
- AMINO ACID and Peptide analysis.
- SUBSTANCE P assay.
- HORMONE and Enzyme assays.
- TISSUE EXTRACTS and Column eluates.
- ELECTROPHORESIS AND CHROMATOGRAPHY is greatly improved by concentration.

You owe it to your lab to order today! Call or write:

**Savant** SAVANT INSTRUMENTS, INC. 221 Park Avenue, Hicksville, NY 11801 • (516) 935-8774